

## Medical Devices Safety Notice

The National Health Regulatory Authority would like to alert all governmental and private healthcare facilities, local agents and distributors that the below medical device:

Device Details	
Device Name	McGRATH™ MAC Video Laryngoscope
Device Model	300-000-000
Serial Number	All serial Numbers prior to 366170*
Manufacturer	Medtronic
Country of Origin	USA
Reference	<a href="#">Link</a>
Reason of Recall	NHRA initiates this FSN due to the risk of battery-related thermal events in McGRATH™ MAC video laryngoscopes, which could lead to battery explosions and potential harm to caregivers or patients.
Action should be taken	Please stop using the above mentioned medical device and contact the authorized representative Cigalah Gulf Medical at <a href="mailto:Maha@cgmed.com">Maha@cgmed.com</a> to take the necessary action for recall.

Your cooperation is highly appreciated in improving health services in the Kingdom of Bahrain.

For more information please contact [Medical\\_Devices@nhra.bh](mailto:Medical_Devices@nhra.bh)